



## Laboratory Single Shaft Slurry Mixing Equipment 150ml For Battery Research

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: MRBEST
- Certification: CE
- Model Number: MR-ZKJB-5
- Minimum Order Quantity: 1 Set
- Price: Negotiable
- Packaging Details: Standard Export Wooden Packing
- Delivery Time: 7-20 working days
- Payment Terms: T/T, Western Union, MoneyGram
- Supply Ability: 500 sets per month



#### Product Specification

- Name: Laboratory Single Shaft Slurry Vacuum Mixing Machine
- Mixing Container: 150ml As Standard, Other Size Is Optional
- Mixing Container Material: Stainless Steel
- Mixing Paddle: Single Shaft
- Vacuum Pump: Built-in Type
- Mixing Speed: 320rpm
- Highlight: **Single Shaft Slurry Mixing Equipment, Laboratory Slurry Mixing Equipment**

for more products please visit us on [lithiumbatterymachine.com](http://lithiumbatterymachine.com)

## Product Description

### Laboratory Single Shaft Slurry Vacuum Mixing Machine 150ml For Battery Research



**Model:** MR-ZKJB-5

Small Lab Vacuum Mixer Machine MR-ZKJB-5 is a desktop vacuum mixer with vacuum mixing with dispersion function. It is widely used for battery slurry mixing, powder and ceramic material mixing under vacuum condition to remove the gas bubble. It is with 150ml mixing container, other volume of mixing container is optional.

#### Features

- Built-in vacuum pump for space saving purpose.
- Volume of mixing container can be optional.
- Mixing time is adjustable.
- Mixing speed is adjustable according to users actual requirement.
- Wide application.

#### Technical Parameter

Model	Small Lab Vacuum Mixer Machine MR-ZKJB-5
Working voltage	AC110V/220V optional
Mixing container	150ml as standard, other size is optional
Mixing container material	Stainless steel
Mixing speed	320rpm
Mixing paddle	Single shaft
Mixing time	0~99.99 hrs settable
Vacuum pump	Built-in type
Vacuum degree	> -80Kpa
Dimension	L330 * W310 * H510mm
Weight	15kg

#### Our Advantages

1. Steady and high quality ! -cause we are a manufacturer of Lithium Battery Equipment, and we only focus on one thing that is make professional Lithium ion Battery Equipment for customer around the world.
2. Friendly use! all MRBEST Battery Equipment with teach manual and support live teaching via media.
3. Reliable factory! My factory located in Shenzhen city of China, support factory visit and video factory visit .
4. 1 year warranty plus lifetime technical support.
5. All the parts use on machine, in stock forever.
6. Powerful technical support, NOT TRADE COMPANY! We are professional on Lithium Battery Equipment, the chief engineer team have more than 20 years experience on lithium battery area, so we can easy do diagnose on machine if customer have question.
7. 100% PATIENCE to all customers! Our patience not only before consult order, also after-sale service, one by one .





+86 19866156608



tina@szmrbest.com



[lithiumbattery-machine.com](http://lithiumbattery-machine.com)

104 Building 6, Second Industrial Zone, Shanmen Community, Yanluo Street, Baoan District, Shenzhen China