



Hand Operated Battery Electrtode Rolling Press Machine For Laboratory

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: MRBEST
- Certification: CE
- Model Number: MR-SG100W
- Minimum Order Quantity: 1 Set
- Price: Negotiable
- Packaging Details: Standard Export Wooden Packing
- Delivery Time: 7-20 working days
- Payment Terms: T/T, Western Union, MoneyGram
- Supply Ability: 500 sets per month



Product Specification

- Product Name: Battery Electrtode Rolling Press Machine
- Roller Diameter: 96mm
- Roller Width: 100mm Other Width Available
- Pressing Thickness: 0~2mm Adjustable
- Rolling Press Accuracy: $\leq \pm 0.005$
- Roller Material: Made From Hardened Tool Steel (HRC > 60)
With Cr Electroplated Coating.
- With finish: > 0.8
- Structure: Heavy Duty Steel Frame, Gears And Handle
- Highlight: **Battery Electrtode Rolling Press Machine,
Laboratory Electrtode Rolling Press Machine**

Product Description

Hand Operated Battery Electrtode Rolling Press Machine for Laboratory



Model: MR-SG100W

MR-SG100W is used for battery electrode rolling press research after coating and drying. It is a bench-top manual rolling press machine with 100mm width roller, and easy to operate.

Features

- Rolling pressing speed can be adjustable.
- High hardness roller, >HRC62
- Forward & backward are available.
- Manual operate, saving power.

Technical Parameter

Model	MR-SG100W
Roller Diameter	φ96
Roller Width	100mm (Other width can be customized)
Roller hardness	>HRC62
Rolling speed	0~10mm/s, adjustable
Weight	45kg

Our Advantages

1. Steady and high quality ! -cause we are a manufacturer of Lithium Battery Equipment, and we only focus on one thing that is make professional Lithium ion Battery Equipment for customer around the world.
2. Friendly use! all MRBEST Battery Equipment with teach manual and support live teaching via media.
3. Reliable factory! My factory located in Shenzhen city of China, support factory visit and video factory visit .
4. 1 year warranty plus lifetime technical support.
5. All the parts use on machine, in stock forever.
6. Powerful technical support, NOT TRADE COMPANY! We are professional on Lithium Battery Equipment, the chief engineer team have more than 20 years experience on lithium battery area, so we can easy do diagnose on machine if customer have question.
7. 100% PATIENCE to all customers! Our patience not only before consult order, also after-sale service, one by one .





+86 19866156608



tina@szmrbest.com



lithiumbattery-machine.com

104 Building 6, Second Industrial Zone, Shanmen Community, Yanluo Street, Baoan District, Shenzhen China