



100mm Width Hot Rolling Press Machine 720W With Winder And Unwinder

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: MRBEST
- Certification: CE
- Model Number: MR-JS100C
- Minimum Order Quantity: 1 Set
- Price: Negotiable
- Packaging Details: Standard Export Wooden Packing
- Delivery Time: 7-20 working days
- Payment Terms: T/T, Western Union, MoneyGram
- Supply Ability: 500 sets per month



Product Specification

- Product Name: Hot Rolling Press Machine With Winder And Unwinder
- Roller: Dia 96 X W100mm, Other Width Available
- Roller Hardness: >HRC62
- Pressing Thickness: 0~2mm Adjustable
- Feeding Speed: Max 10mm/s Adjustable
- Heating Temperature: RT~ 120 Adjustable digital Display
- Working Voltage: AC110/220V Optional Upon Request
- Power: 720W
- Highlight: **100mm Hot Rolling Press Machine,
Hot Rolling Press Machine 720W**

Product Description

100mm Width Hot Rolling Press Machine With Winder And Unwinder



Model: MR-JS100C

Applications:

MR-JS100C hot rolling press machine with Winder and Unwinder is used for lab lithium battery electrode calendaring after mental foils coating and drying process. it can equipped with winding & unwinding device for battery electrode after calendaring. Two sets of digital dial indicator for pressing thickness setting easily, heating temp is adjustable for RT~150 .

Main features

- It is equipped with winding and unwinding device for battery electrode.
- Built-in type heating tube is adopted designing inside of the rollers for heating effect more evenly.
- Adjustable pressing thickness, digital dial indicators with memory function when shutdown,
- Digital display speed controller, easy to control and adjust the pressing speed.
- Vertical rollers installation design for convenient feeding and pressing, forward *backward are both available.
- Heating temperature is adjustable and can be controlled separately by thermostats.

Specifications:

Model	MR-JS100C
Voltage	AC110/220V optional
Power	720W
Roller size	Φ100 * L100mm
Roller hardness	>HRC62
Pressing thickness	0.01~2mm adjustable (0.001mm optional)
Heating temp	RT~150 adjustable, digital display
Pressing speed	0~10mm/s, digital display
Dimension	L670 * W225 * H300mm
Weight	63kg

Our Advantages

1. Steady and high quality ! -cause we are a manufacturer of Lithium Battery Equipment, and we only focus on one thing that is make professional Lithium ion Battery Equipment for customer around the world.
2. Friendly use! all MRBEST Battery Equipment with teach manual and support live teaching via media.
3. Reliable factory! My factory located in Shenzhen city of China, support factory visit and video factory visit .
4. 1 year warranty plus lifetime technical support.
5. All the parts use on machine, in stock forever.
6. Powerful technical support, NOT TRADE COMPANY! We are professional on Lithium Battery Equipment, the chief engineer team have more than 20 years experience on lithium battery area, so we can easy do diagnose on machine if customer have question.
7. 100% PATIENCE to all customers! Our patience not only before consult order, also after-sale service, one by one .



Shenzhen Meirui Zhida Technology Co., Ltd.



+86 19866156608



tina@szmrbest.com



lithiumbattery-machine.com

104 Building 6, Second Industrial Zone, Shanmen Community, Yanluo Street, Baoan District, Shenzhen China