



MR-CP60 Manual Lithium Battery Making Machine Coin Cell Disc Cutter

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: MRBEST
- Certification: CE
- Model Number: MR-CP60
- Minimum Order Quantity: 1 Set
- Price: Negotiable
- Packaging Details: Standard Export Wooden Packing
- Delivery Time: 7-20 working days
- Payment Terms: T/T, Western Union, MoneyGram
- Supply Ability: 500 sets per month



Product Specification

- Name: Coin Cell Disc Cutter
- Cutting Die: Optional 3 - 24mm Customized Dies Are Available Upon Request
- Structure: Adopt Aluminum Alloy Corrosion-resistant Material And Never Rusts.
- Max. Stroke Length: 25mm
- Product Dimensions: 110mm(L) X 160mm(W) X 350mm(H)
- Net Weight: 6kg
- Highlight: **Manual Lithium Battery Making Machine, Coin Cell Lithium Battery Making Machine, Manual Coin Cell Disc Cutter**

Product Description

Battery Making Machine Coin Cell Disc Cutter For Battery Lab Research



Model: MR-CP60

The Coin Cell Disc Cutter Machine is used for coin cell electrode punching with customized cutting die. It's operation is simple with light weight, The resulting discs can be used as electrode/separator discs of split test cell and coin cell for battery researches.

Cutting Die: 3 - 24mm Optional

Specifications:

Machine Model	Coin Cell Disc Cutter MR-CP60
Operating	Manual
Punching size	Optional 3 - 24mm customized dies are available upon request
Application Notes	For make CR2016, CR2025 or/and CR2032 coin cell batteries, we will need at least three cutting dies: A 15mm B 16mm C 19 mm (for separator to fully separate the electrodes) For CR2450 coin cell: 20mm cutting die can be used for cutting electrode discs and the 23mm cutting die for separator disc.
Punching stroke	25mm
Punching thickness	0.005-0.5mm electrode and separator
Advantages	Simple operation Light weight Can both cut electrode and separator
Dimension	110*150*235mm
Weight	About 6kg
Warranty	One Year limited warranty with lifetime support

Our Advantages

1. Steady and high quality ! -cause we are a manufacturer of Lithium Battery Equipment, and we only focus on one thing that is make professional Lithium ion Battery Equipment for customer around the world.
2. Friendly use! all MRBEST Battery Equipment with teach manual and support live teaching via media.
3. Reliable factory! My factory located in Shenzhen city of China, support factory visit and video factory visit .
4. 1 year warranty plus lifetime technical support.
5. All the parts use on machine, in stock forever.
6. Powerful technical support, NOT TRADE COMPANY! We are professional on Lithium Battery Equipment, the chief engineer team have more than 20 years experience on lithium battery area, so we can easy do diagnose on machine if customer have question.
7. 100% PATIENCE to all customers! Our patience not only before consult order, also after-sale service, one by one .



Shenzhen Meirui Zhida Technology Co., Ltd.



+86 19866156608



tina@szmrbest.com



lithiumbattery-machine.com

104 Building 6, Second Industrial Zone, Shanmen Community, Yanluo Street, Baoan District, Shenzhen China