



Laboratory Vacuum Oven

Our Product Introduction

for more products please visit us on lithiumbatterymachine.com

Basic Information

- Place of Origin: China
- Brand Name: MRBEST
- Certification: CE
- Model Number: MR-DZF-6050
- Minimum Order Quantity: 1 Set
- Price: Negotiable
- Packaging Details: Standard Export Wooden Packing
- Delivery Time: 7-20 working days
- Payment Terms: T/T, Western Union, MoneyGram
- Supply Ability: 500 sets per month



Product Specification

- Product Name: Laboratory Vacuum Oven
- Temperature: RT+10 ~ 200°C Is Standard, RT+10 ~ 250°C Is Optional
- Volume: 50L (other Volume Optional)
- Input Voltage: 1400W
- Vacuum Degree: 133Pa
- Dimensions: Oven: 710 X 545 X 542 Mm
- Temperature Controller: Built-in Over-temperature Protection
- Temperature Control Accuracy: +/-1°C
- Highlight: **laboratory lithium battery making machine,
laboratory lithium ion battery making machine**

Laboratory Vacuum Oven For Battery Material/Electrode Baking



Model: MR-DZF-6050

1. The DZF-6050 Vacuum Drying Oven is mainly used For scientific research units, specialty universities, laboratories, industrial mining enterprises etc and production field for materials drying and heat treatment under the vacuum state.
2. The DZF-6050 Vacuum Drying Oven for lithium battery field mainly used for lithium battery materials, electrode and battery core vacuum drying.
3. Use DZF-6050 Vacuum drying oven for material heating drying under vacuum states has advantages :
 - 3.1 Can reduce the drying temperature
 - 3.2 Avoid some objects heating to oxidation
 - 3.3 Avoid heating air to kill living cells
 - 3.4 No dust damage

Specifications:

Model	Laboratory Drying Oven MR-DZF-6050
Source voltage	110VAC or 220VAC +/-10% selectable
Power	1400W
Warranty	One Year Limited Manufacture Warranty with Life Time Support
Temperature range	RT+10-250

Temperature fluctuation	±0.5
Vacuum degree	-0.1MPa
Inner box size(L*W*H)	415*345*370mm
Dimension(L*W*H)	710*560*550mm
Package size	650*850*700mm
Temperature controller	1 High precision microcomputer temperature controller, accurate temperature control 2 Working temperature can be set and working time can be set for 0-9999min. 3 When time reach, sound alarm, and when any error, also alarm. function
Inner box material	304 stainless steel
Outside shell material	High quality cold-rolled steel plate electrostatic spray
Box door	Toughened Bulletproof double glass door
Safety device	With over temperature protection, grounding etc
Shelf	2pcs 1 Solid stainless steel shelf, guarantee the oven temperature 2 Mesh stainless steel shelf can be also chosen

Our Advantages

1. Steady and high quality ! -cause we are a manufacturer of Lithium Battery Equipment, and we only focus on one thing that is make professional Lithium ion Battery Equipment for customer around the world.
2. Friendly use! all MRBEST Battery Equipment with teach manual and support live teaching via media.
3. Reliable factory! My factory located in Shenzhen city of China, support factory visit and video factory visit .
4. 1 year warranty plus lifetime technical support.
- 5.All the parts use on machine, in stock forever.
6. Powerful technical support, NOT TRADE COMPANY! We are professional on Lithium Battery Equipment, the chief engineer team have more than 20 years experience on lithium battery area, so we can easy do diagnose on machine if customer have question.
7. 100% PATIENCE to all customers! Our patience not only before consult order, also after-sale service, one by one .



Shenzhen Meirui Zhida Technology Co., Ltd.



+86 19866156608



tina@szmrbest.com



lithiumbatteryemachine.com

104 Building 6, Second Industrial Zone, Shanmen Community, Yanluo Street, Baoan District, Shenzhen China