



Lithium Ion Battery Manufacturing Machine Coin Cell Battery Equipment

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: MRBEST
- Certification: CE
- Model Number: MRBEST-Coin Cell
- Minimum Order Quantity: 1 Set
- Price: Negotiable
- Packaging Details: Standard Export Wooden Packing
- Delivery Time: 7-20 working days
- Payment Terms: T/T, Western Union, MoneyGram
- Supply Ability: 500 sets per month



Product Specification

- Product Name: Coin Cell Battery Equipment
- Application: For Battery Research
- Suitable Battery: Coin Battery
- Size: Accepted Customized Size
- Manufacturer: YES
- Direct Factory: With Competitive Price
- High Light: Battery Manufacturing Machine, Battery Making Machine Price, Lithium Battery Making Machine
- Highlight: **Lithium ion Battery Manufacturing Machine, Coin Cell Battery Manufacturing Machine, Lithium Ion Coin Cell Battery Equipment**

Coin Cell Battery Equipment

Coin Cell Preparation



Shenzhen Meirui Zhida Technology Co., Ltd
Tel/Fax: +86 755 23696396 meiruitop@gmail.com www.mrbest.en.alibaba.com

Coin Cell Battery alternatively known as a button cell or watch battery, a coin cell battery is a small, circular, stainless steel disc that provides power to portable devices, like watches and computers.

How to Make a Coin Cell?

MRBEST can provide the complete line machines to make a coin cell, followings are the list of Coin Cell Battery Equipment:

No.	Coin Cell Battery Equipment
1	<u>Tube Furnace/Muffle Furnace</u>
2	Planetary Ball Mill
3	Vacuum Mixer
4	Film Coater Machine
5	Roll Press Machine
6	Coin Cell Punching Machine
7	Vacuum Drying Oven
8	Glove Box Machine
9	Coin Cell Crimper
10	Battery Tester Machine
11	E-mail: tina@szmrbest.com

Our Advantages

1. Steady and high quality ! -cause we are a manufacturer of Lithium Battery Equipment, and we only focus on one thing that is make professional Lithium ion Battery Equipment for customer around the world.
2. Friendly use! all MRBEST Battery Equipment with teach manual and support live teaching via media.
3. Reliable factory! My factory located in Shenzhen city of China, support factory visit and video factory visit .
4. 1 year warranty plus lifetime technical support.

5.All the parts use on machine, in stock forever.

6. Powerful technical support, NOT TRADE COMPANY! We are professional on Lithium Battery Equipment, the chief engineer team have

more than 20 years experience on lithium battery area, so we can easy do diagnose on machine if customer have question.

7. 100% PATIENCE to all customers! Our patience not only before consult order, also after-sale service, one by one .



Shenzhen Meirui Zhida Technology Co., Ltd.



+86 19866156608



tina@szmrbest.com



lithiumbattery-machine.com

104 Building 6, Second Industrial Zone, Shanmen Community, Yanluo Street, Baoan District, Shenzhen China