



Vacuum Battery Slurry Mixer Max 1400RPM Speed For Laboratory Research

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: MRBEST
- Certification: CE
- Model Number: MR-XJB-4
- Minimum Order Quantity: 1 Set
- Price: Negotiable
- Packaging Details: Standard Export Wooden Packing
- Delivery Time: 7-20 working days
- Payment Terms: T/T, Western Union, MoneyGram
- Supply Ability: 500 sets per month



Product Specification

- Name: Battery Slurry Vacuum Mixer For Laboratory Research
- Mixing Container: 500ml As Standard, 150ml & 250ml Is Optional
- Mixing Container Material: Stainless Steel
- Mixing Speed: Max 1400RPM
- Mixing Paddle: Single Shaft
- Vacuum Pump: Built-in
- Dimension: L330 * W320 * H570mm
- Weight: 35KGS
- Highlight: **Vacuum Battery Slurry Mixer,
Laboratory Battery Slurry Mixer**

for more products please visit us on lithiumbatteryachine.com

Product Description

Battery Slurry Vacuum Mixer For Laboratory Research



Model: MR-XJB-4

Product Introduction

MR-XJB-4 is a desktop vacuum mixer with vacuum mixing and dispersion function. It is widely used for battery slurry mixing, powder and ceramic material mixing under vacuum condition to remove the gas bubble. It is with 150ml mixing container, other volume of mixing container is optional.

Features

- Built-in vacuum pump for space saving purpose.
- Volume of mixing container can be optional.
- Mixing time is adjustable.
- Mixing speed is adjustable according to users actual requirement.
- Wide application.

Specifications

Name	Battery Slurry Vacuum Mixer For Laboratory Research
Model	MR-XJB-4
Power	AC110V/220V optional upon request
Mixing container	500ml as standard, other size is optional
Mixing container material	Stainless steel
Mixing speed	Max 1400rpm, 3 segments programmable
Mixing paddle	Single shaft
Mixing time	0~99.99 hrs settable
Vacuum pump	Built-in type
Vacuum degree	-0.08Mpa
Dimension	L330 * W320 * H570mm
Weight	35kg

Standard Package

Item No.	Description	Quantity
MR-XJB-4	Compact Dual-Shaft Planetary Vacuum Mixer	1

Stainless Steel Vacuum Sealed Container	1 stainless steel vacuum sealed container are included (optional 150ml, 250ml or 500ml).	1
---	--	---

Our Advantages

1. Steady and high quality ! -cause we are a manufacturer of Lithium Battery Equipment, and we only focus on one thing that is make professional Lithium ion Battery Equipment for customer around the world.
2. Friendly use! all MRBEST Battery Equipment with teach manual and support live teaching via media.
3. Reliable factory! My factory located in Shenzhen city of China, support factory visit and video factory visit .
4. 1 year warranty plus lifetime technical support.
- 5.All the parts use on machine, in stock forever.
6. Powerful technical support, NOT TRADE COMPANY! We are professional on Lithium Battery Equipment, the chief engineer team have more than 20 years experience on lithium battery area, so we can easy do diagnose on machine if customer have question.
7. 100% PATIENCE to all customers! Our patience not only before consult order, also after-sale service, one by one .



Shenzhen Meirui Zhida Technology Co., Ltd.



+86 19866156608



tina@szmrbest.com



lithiumbatteryachine.com

104 Building 6, Second Industrial Zone, Shanmen Community, Yanluo Street, Baoan District, Shenzhen China