



Compact Coin Cell Battery Electrode Disc Cutter For Precise Round Disc Cutting

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: MRBEST
- Certification: CE
- Model Number: MR-CP60
- Minimum Order Quantity: 1 Set
- Price: Negotiable
- Packaging Details: Standard Export Wooden Packing
- Delivery Time: 7-20 working days
- Payment Terms: T/T, Western Union, MoneyGram
- Supply Ability: 500 sets per month



Product Specification

- Name: Coin Cell Electrode Disc Cutter
- Cutting Die: Optional 3 - 24mm Customized Dies Are Available Upon Request
- Structure: Adopt Aluminum Alloy Corrosion-resistant Material And Never Rusts.
- Max. Stroke Length: 25mm
- Product Dimensions: 110mm(L) X 160mm(W) X 350mm(H)
- Net Weight: 6kg
- Highlight: **Compact Battery Electrode Disc Cutter ,
Precise Battery Electrode Disc Cutter ,
Compact Coin Cell Electrode Disc Cutter**

Product Description

Lithium Ion Battery Assembly Line Coin Cell Electrode Disc Cutter



Model: MR-CP60

This compact disc cutter is used for precise round disc cutting. Standard cutting die $\phi 3\sim 24\text{mm}$ can be optional. It is available to cut various metal and foil thickness $0.005 \sim 0.5\text{mm}$, especially it's a perfect tool for cutting coin cell electrode and separator. It is small size can be put into glovebox through antechamber.

Cutting Die: 3 - 24mm Optional

Main features

- Different size of cutting die can be optional. Standard die $\phi 3\sim 24\text{mm}$.
- Fast to change the cutting die.
- High accuracy, no burrs and dents.
- Various of battery material with thickness $0.005\sim 0.5\text{mm}$ can be cut.
- Small dimension, convenient operation, quick and flexible use into glovebox.
- The machine body is made of anti-corrosion material that is anti-rust.

Specifications:

Cutting die	$\phi 3\sim 24\text{mm}$ optional, customized size is available
Cutting material thickness	$0.005\sim 0.5\text{mm}$ various of battery material
Punching stroke	25mm
Operation mode	Manual
Dimension	L110xW150xH220mm
Weight	10kg

Our Advantages

1. Steady and high quality ! -cause we are a manufacturer of Lithium Battery Equipment, and we only focus on one thing that is make professional Lithium ion Battery Equipment for customer around the world.
2. Friendly use! all MRBEST Battery Equipment with teach manual and support live teaching via media.
3. Reliable factory! My factory located in Shenzhen city of China, support factory visit and video factory visit .
4. 1 year warranty plus lifetime technical support.
5. All the parts use on machine, in stock forever.
6. Powerful technical support, NOT TRADE COMPANY! We are professional on Lithium Battery Equipment, the chief engineer team have

more than 20 years experience on lithium battery area, so we can easy do diagnose on machine if customer have question.
7. 100% PATIENCE to all customers! Our patience not only before consult order, also after-sale service, one by one .



Shenzhen Meirui Zhida Technology Co., Ltd.



+86 19866156608



tina@szmrbest.com



lithiumbatterymachine.com

104 Building 6, Second Industrial Zone, Shanmen Community, Yanluo Street, Baoan District, Shenzhen China