



High Precision Battery Manufacturing Machine Electric Sealer Coin Cell Crimper Machine

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: MRBEST
- Certification: CE
- Model Number: MR-DF160
- Minimum Order Quantity: 1 Set
- Price: Negotiable
- Packaging Details: Standard Export Wooden Packing
- Delivery Time: 7-20 working days
- Payment Terms: T/T, Western Union, MoneyGram
- Supply Ability: 500 sets per month



Product Specification

- Name: Electric Sealer Coin Cell Crimper Machine
- Application: For Lab Button Battery Sealing
- Standard Crimping Die: CR20xx Crimping Die
- Operation Mode: Electric
- Size: L385xW160xH325mm
- Weight: 25KGS
- Highlight: **High Precision Battery Manufacturing Machine,
Electric Battery Manufacturing Machine,
Electric Coin Cell Crimper Machine**

for more products please visit us on lithiumbatteryemachine.com

Product Description

Electric Sealer Coin Cell Crimper Machine for Crimping CR20XX Coin Cells



Model: MR-DF160

Profile:

It is a compact electric crimper for various of coin cell crimping, such as CR2016, CR2025, CR2032 and etc. CR20xx series crimping die is standard, the other crimping die set can be optional. It is touch screen controlled with high automation. Compact size makes it easy to enter into glovebox. It is a perfect device for coin cell research for universities, research institutes, factories and etc.

Features :

- Stainless steel crimping die, anti-corrosion and durable.
- Touch screen operated, high automation.
- Compact size, easy to use inside of glovebox.
- Various of crimping die is optional.
- High precision mould design to ensure stable sealing.
- It's designed with crimping height limitation to avoid damage.

Technical Specifications:

Brand name	MRBEST
Power supply	AC110V/220V optional
Standard crimping die	CR20xx crimping die, other size crimping die is optional
Coin cell disassemble die	Optional
Operation	Electric mode
Dimension	L385xW160xH325mm
Weight	32kg

Our Advantages

1. Steady and high quality ! -cause we are a manufacturer of Lithium Battery Equipment, and we only focus on one thing that is make professional Lithium ion Battery Equipment for customer around the world.
2. Friendly use! all MRBEST Battery Equipment with teach manual and support live teaching via media.
3. Reliable factory! My factory located in Shenzhen city of China, support factory visit and video factory visit .
4. 1 year warranty plus lifetime technical support.
5. All the parts use on machine, in stock forever.
6. Powerful technical support, NOT TRADE COMPANY! We are professional on Lithium Battery Equipment, the chief engineer team have more than 20 years experience on lithium battery area, so we can easy do diagnose on machine if customer have question.
7. 100% PATIENCE to all customers! Our patience not only before consult order, also after-sale service, one by one .



Shenzhen Meirui Zhida Technology Co., Ltd.

