

5V 10mA / 20mA / 50mA Battery Testing Equipment 8 Channels BTS-4000 Series

Our Product Introduction

for more products please visit us on lithiumbatteryemachine.com

Basic Information

- Place of Origin: China
- Brand Name: MRBEST
- Certification: CE
- Model Number: CT-4008T-5V10mA/5V20mA/ 5V50mA
- Minimum Order Quantity: 1 Set
- Price: Negotiable
- Packaging Details: Standard Export Wooden Packing
- Delivery Time: 7-20 working days
- Payment Terms: T/T, Western Union, MoneyGram
- Supply Ability: 500 sets per month



Product Specification

- Name: BTS-4000 Series 5V10mA 5V20mA 5V50mA Battery Tester
- Per Channel Current Range: Range 1:5 μ A ~ 1 MA; Range 2: 1mA ~ 25mA; Range 3: 25 MA~50mA
- Channel Features: Constant Current Source And Constant Voltage Source With Independent Pairs Of Closed-loop Structure
- Channel Control Mode: Independent Control
- Per Channel Voltage Range: 25mV~5V
- Single-channel Maximum Output Power: 0.25W
- Highlight: 5V 50mA Battery Testing Equipment, 5V 20mA Battery Testing Equipment

Product Description

BTS-4000 Series 5V10mA/ 5V20mA/ 5V50mA Battery Tester**Model:** CT-4008T-5V10mA/ 5V20mA/ 5V50mA

BTS-4000 series entered into market in the year of 2008. Now it is one of the most successful and most popular battery testing system in China and World. Thousands of battery manufacturers, battery research institutes and laboratories use it as their primary testing partner for their various testing activities.

Due to the rapid data acquisition frequency, BTS-4000 could be applied to super capacitor or EDLC (Electronic Double Layer Capacitor) test. Through our advanced software, the capacitance, DCIR(Direct Current Internal Resistance) could be easily calculated.

Smaller, More Powerful

Higher Accuracy. More Cost-effective. We believe it can improve your battery test.

Triple Range

Multi range improve your test accuracy! Automatic switched according to the input current value.

CT-4008T-5V10mA-164

Range1: 5 μ A~1mA; Range2: 1mA~5mA; Range3: 5mA~10mA; Accuracy: 0.05% FS

CT-4008T-5V20mA-164


Range1: 5 μ A~1mA; Range2: 1mA~10mA; Range3: 10mA~20mA; Accuracy: 0.05% FS

CT-4008T-5V50mA-164

Range1: 5 μ A~1mA; Range2: 1mA~25mA; Range3: 25mA~50mA; Accuracy: 0.05% FS

High Accuracy Battery Cycler

Eight-channels battery charge and discharge testing system




Multi Channels can be Supported!


Neware design the practical rack that customer can put the tester in this rack at the best location. (Maximum 20 set of mA Neware battery testers can be put in one rack)

Different Batteries can be test


- * Coin Cell
- * Three-electrode System
- * Flow Battery
- * Lithium-Air Batteries
- * Sodium-ion Battery
- * Li--Sulfur Battery
- * Solid-State Battery
- * Other types Batteries




Kinds of Testing Steps




Different Expression Supported



Customized Parameter Exported



Multi Range Current Switched



Amazing Testing Cable



Beautiful Testing Clamps



Coin Cell Specialty



Shenzhen Meirui Zhida Technology Co., Ltd.



+86 19866156608



tina@szmrbest.com



lithiumbatterymachine.com

104 Building 6, Second Industrial Zone, Shanmen Community, Yanluo Street, Baoan District, Shenzhen China